

# ZNF256 Antibody



PACO13322

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZNF256

**Uniprot**

Q9Y2P7

**Synonyms:**

55kd;

**Immunogen:**

Human ZNF256.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A