

ZNF35 Antibody



PACO13327

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF35

Uniprot

P13682

Synonyms:

zinc finger protein 35; ZNF35; HF.10; HF10; Zfp105 ;

Immunogen:

Human ZNF35.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A