

ZNF385D Antibody



PACO13328

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF385D

Uniprot

Q9H6B1

Synonyms:

zinc finger protein 385D; ZNF385D; FLJ22419; ZNF659 ;

Immunogen:

Human ZNF385D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A