

ZNF443 Antibody



PACO13331

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ZNF443

Uniprot

Q9Y2A4

Synonyms:

zinc finger protein 443; ZNF443; ZK1 ;

Immunogen:

Human ZNF443.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A