

# ZNF622 Antibody



PACO13337

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

ZNF622

**Uniprot**

Q969S3

**Synonyms:**

zinc finger protein 622; ZNF622; MGC17552; MGC2485; ZPR9 ;

**Immunogen:**

Human ZNF622.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A