

PACO13345

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ZXDC

Uniprot

Q2QGD7

Synonyms:

ZXD family zinc finger C; ZXDC; DKFZp547N024; FLJ13861; FLJ36477; FLJ37830; MGC11349; ZIDL ;

Immunogen:

Human ZXDC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A