

ACTA1 Antibody



PACO13355

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ACTA1

Uniprot

P68133

Synonyms:

actin, alpha 1, skeletal muscle; ACTA1; ACTA; ASMA; CFTD; CFTD1; CFTDM; MPFD; NEM1; NEM2; NEM3 ; ; Alpha-actin-1

Immunogen:

Human ACTA1-skeletal-Muscle-Actin-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A