

CD80 Antibody



PACO13375

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD80

Uniprot

P33681

Synonyms:

CD80 molecule; CD80; CD28LG; CD28LG1; LAB7 ;

Immunogen:

Human B7-1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A