

# ATG14 Antibody



PACO13377

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ATG14

**Uniprot**

Q6ZNE5

**Synonyms:**

KIAA0831; KIAA0831; ATG14; MGC126291; MGC126292 ;

**Immunogen:**

Human Barkor-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A