

MUC1 Antibody



PACO13388

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MUC1

Uniprot

P15941

Synonyms:

mucin 1, cell surface associated; MUC1; CD227; EMA; H23AG; MAM6; PEM; PEMT; PUM ; MUC1

Immunogen:

Human CA15-3, MUC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A