

EIF4E Antibody



PACO13389

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

EIF4E

Uniprot

P06730

Synonyms:

eukaryotic translation initiation factor 4E; EIF4E; CBP; EIF4E1; EIF4EL1; EIF4F; MGC111573 ;

Immunogen:

Human CBP, EIF4E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A