

PACO13393

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CCR2

Uniprot

P41597

Synonyms:

chemokine (C-C motif) receptor 2; CCR2; CC-CKR-2; CCR2A; CCR2B; CD192; CKR2; CKR2A; CKR2B; CMKBR2; FLJ78302; MCP-1-R; MGC103828; MGC111760; MGC168006 ;

Immunogen:

Human CCR2a-specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A