

CCT6A Antibody



PACO13394

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CCT6A

Uniprot

P40227

Synonyms:

chaperonin containing TCP1, subunit 6A (zeta 1); CCT6A; CCT-zeta; CCT-zeta-1; CCT6; Cctz; HTR3; MGC126214; MGC126215; MoDP-2; TCP-1-zeta; TCP20; TCPZ; TTCP20 ;

Immunogen:

Human CCT6A-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A