

# CDKN1C Antibody



PACO13400

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDKN1C

**Uniprot**

P49918

**Synonyms:**

cyclin-dependent kinase inhibitor 1C (p57, Kip2); CDKN1C; BWCR; BWS; KIP2; WBS; p57 ;

**Immunogen:**

Human CDKN1C-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A