

# CLDN4 Antibody



PACO13404

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CLDN4

**Uniprot**

O14493

**Synonyms:**

claudin 4; CLDN4; CPE-R; CPER; CPETR; CPETR1; WBSCR8; hCPE-R;

**Immunogen:**

Human CLDN4-specific.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A