

# COL16A1 Antibody



PACO13406

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COL16A1

**Uniprot**

Q07092

**Synonyms:**

collagen, type XVI, alpha 1; COL16A1;447AA; FP1572 ;

**Immunogen:**

Human COL16A1-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A