

GJB2 Antibody



PACO13407

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GJB2

Uniprot

P29033

Synonyms:

gap junction protein, beta 2, 26kDa; GJB2; CX26; DFNA3; DFNA3A; DFNB1; DFNB1A; HID; KID; NSRD1; PPK ; Connexin 26

Immunogen:

Human Connexin-26.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A