

COX4I1 Antibody



PACO13409

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

COX4I1

Uniprot

P13073

Synonyms:

cytochrome c oxidase subunit IV isoform 1; COX4I1; COX4; COXIV; MGC72016 ;

Immunogen:

Human COXIV.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A