

CPXM2 Antibody



PACO13411

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CPXM2

Uniprot

Q8N436

Synonyms:

carboxypeptidase X (M14 family), member 2; CPXM2; UNQ676 ;

Immunogen:

Human CPXM2-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A