

CREB3L1 Antibody



PACO13412

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CREB3L1

Uniprot

Q96BA8

Synonyms:

cAMP responsive element binding protein 3-like 1; CREB3L1; OASIS ;

Immunogen:

Human CREB3L1, OASIS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A