

CSTL1 Antibody



PACO13418

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CSTL1

Uniprot

Q9H114

Synonyms:

cystatin-like 1; CSTL1; RCET11; dJ322G13.4 ;

Immunogen:

Human CSTL1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A