

MARK3 Antibody



PACO13419

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MARK3

Uniprot

P27448

Synonyms:

MAP/microtubule affinity-regulating kinase 3; MARK3; CTAK1; KP78; PAR1A ;

Immunogen:

Human CTAK1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A