

CYP4Z1 Antibody



PACO13423

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CYP4Z1

Uniprot

Q86W10

Synonyms:

cytochrome P450, family 4, subfamily Z, polypeptide 1; CYP4Z1; CYP4A20 ;

Immunogen:

Human CYP4Z1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A