

DLL1 Antibody



PACO13434

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DLL1

Uniprot

O00548

Synonyms:

delta-like 1 (Drosophila); DLL1; DELTA1; DL1; Delta ;

Immunogen:

Human DLL1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A