

# DCTN2 Antibody



PACO13440

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DCTN2

**Uniprot**

A8K8J9

**Synonyms:**

dynactin 2 (p50); DCTN2; DCTN50; DYNAMITIN; RBP50 ;

**Immunogen:**

Human dynactin-2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A