

EFNA4 Antibody



PACO13442

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EFNA4

Uniprot

P52798

Synonyms:

ephrin-A4; EFNA4; EFL4; EPLG4; LERK4; MGC125826 ; Ephrin A4; LERK-4

Immunogen:

Human EFNA4-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A