

SLC29A1 Antibody



PACO13444

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC29A1

Uniprot

Q99808

Synonyms:

solute carrier family 29 (nucleoside transporters), member 1; SLC29A1; ENT1; MGC1465; MGC3778 ;

Immunogen:

Human ENT1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A