

EPHA1 Antibody



PACO13445

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EPHA1

Uniprot

P21709

Synonyms:

EPH receptor A1; EPHA1; EPH; EPHT; EPHT1; MGC163163 ;

Immunogen:

Human EPHA1-special.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A