

# MAPK3 Antibody



PACO13446

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat, Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAPK3

**Uniprot**

P27361

**Synonyms:**

mitogen-activated protein kinase 3; MAPK3; ERK1; HS44KDAP; HUMKER1A; MGC20180; P44ERK1; P44MAPK; PRKM3 ; ERK1/2; MAPK3

**Immunogen:**

Human ERK1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A