

EMR1 Antibody



PACO13450

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ADGRE1

Uniprot

Q14246

Synonyms:

EGF-like module containing, mucin-like, hormone receptor-like sequence 1; Emr1; DD7A5-7; EGF-TM7; F4/80; Gpf480; Ly71; TM7LN3 ;

Immunogen:

Human F4-80.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A