

FCGRT Antibody



PACO13452

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FCGRT

Uniprot

P55899

Synonyms:

Fc fragment of IgG, receptor, transporter, alpha; FCGRT; FCRN; alpha-chain ;

Immunogen:

Human FCGRT-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A