

FETUB Antibody



PACO13453

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FETUB

Uniprot

Q9UGM5

Synonyms:

fetuin B; FETUB;16G2; Gugu; IRL685 ;

Immunogen:

Human fetuin-B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A