

LGALS1 Antibody



PACO13457

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LGALS1

Uniprot

P09382

Synonyms:

lectin, galactoside-binding, soluble, 1; LGALS1; DKFZp686E23103; GAL1; GBP ; Galectin-1, S-Lac lectin 1, HPL, HBL, HLBP14

Immunogen:

Human GAL1, LGALS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A