

PACO13461

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GBA

Uniprot

P04062

Synonyms:

glucosidase, beta; acid, pseudogene; GBAP; MGC104662; MGC40022; MGC8614 ;

Immunogen:

Human GBAP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A