

HBM Antibody



PACO13476

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HBM

Uniprot

Q6B0K9

Synonyms:

hemoglobin, mu; HBM; HBAP2 ;

Immunogen:

Human HBM-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A