

HDAC5 Antibody



PACO13477

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HDAC5

Uniprot

Q9UQL6

Synonyms:

histone deacetylase 5; HDAC5; FLJ90614; HD5; NY-CO-9 ;

Immunogen:

Human HDAC5-specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A