

# HSP90AA1 Antibody



PACO13484

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

HSP90AA1

**Uniprot**

P07900

**Synonyms:**

heat shock protein 90kDa alpha (cytosolic), class A member 1; HSP90AA1; FLJ31884; HSP86; HSP89A; HSP90A; HSP90N; HSPC1; HSPCA; HSPCAL1; HSPCAL4; HSPN; Hsp89; Hsp90; LAP2 ;

**Immunogen:**

Human HSP90.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A