

# ING1 Antibody



PACO13491

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ING1

**Uniprot**

Q9UK53

**Synonyms:**

inhibitor of growth family, member 1; ING1; p24ING1c; p33; p33ING1; p33ING1b; p47; p47ING1a ;

**Immunogen:**

Human ING1-specific.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A