

DYNC1H1 Antibody



PACO13500

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

DYNC1H1

Uniprot

Q14204

Synonyms:

dynein, cytoplasmic 1, heavy chain 1; DYNC1H1; DHC1; DHC1a; DKFZp686P2245; DNCH1; DNCL; DNECL; DYHC; Dnchc1; HL-3; KIAA0325; p22 ;

Immunogen:

Human KIAA0325.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A