

# KLK3 Antibody



PACO13501

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KLK3

**Uniprot**

P07288

**Synonyms:**

kallikrein-related peptidase 3; KLK3; APS; KLK2A1; PSA; hK3 ; Prostate-specific antigen; KLK3

**Immunogen:**

Human KLK3, PSA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A