

LAMC3 Antibody



PACO13504

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LAMC3

Uniprot

Q9Y6N6

Synonyms:

laminin, gamma 3; LAMC3; DKFZp434E202 ;

Immunogen:

Human LAMC3-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A