

PACO13506

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RPSA

Uniprot

P08865

Synonyms:

ribosomal protein SA; RPSA;37LRP;67LR; LAMBR; LAMR1; LRP; p40 ; Laminin receptor; 37/67 kDa laminin receptor; LRP/LR; LBP/p40

Immunogen:

Human LAMR1, RPSA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A