

MAP2 Antibody



PACO13514

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MAP2

Uniprot

P11137

Synonyms:

microtubule-associated protein 2; MAP2; DKFZp686l2148; MAP2A; MAP2B; MAP2C ;

Immunogen:

Human MAP2B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A