

MARCH5 Antibody



PACO13515

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

38412

Uniprot

Q9NX47

Synonyms:

membrane-associated ring finger (C3HC4) 5; MARCH5; FLJ20445; MARCH-V; MITOL; RNF153 ;

Immunogen:

Human MARCH5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A