

MARCH8 Antibody



PACO13516

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

39508

Uniprot

Q5T0T0

Synonyms:

membrane-associated ring finger (C3HC4) 8; MARCH8; MARCH-VIII; MIR; RNF178; c-MIR ;

Immunogen:

Human MARCH8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A