

MEF2B Antibody



PACO13518

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MEF2B

Uniprot

Q02080

Synonyms:

LOC729991-MEF2B readthrough transcript; LOC729991-MEF2B; MEF2B ;

Immunogen:

Human MEF2B-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A