

PACO13521

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MME

Uniprot

P08473

Synonyms:

membrane metallo-endopeptidase; MME; CALLA; CD10; DKFZp686O16152; MGC126681; MGC126707; NEP ; Nepriylsin; Enkephalinase; Atriopeptidase

Immunogen:

Human MME, CD10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A