

MYOM1 Antibody



PACO13525

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYOM1

Uniprot

P52179

Synonyms:

myomesin 1, 185kDa; MYOM1; MGC134946; MGC134947; SKELEMIN ;

Immunogen:

Human MYOM1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A