

# NARFL Antibody



PACO13531

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NARFL

**Uniprot**

Q9H6Q4

**Synonyms:**

nuclear prelamin A recognition factor-like; NARFL; FLJ21988; HPRN; IOP1; LET1L; PRN ;

**Immunogen:**

Human NARFL-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A