

CDH2 Antibody



PACO13532

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CDH2

Uniprot

P19022

Synonyms:

cadherin 2, type 1, N-cadherin (neuronal); CDH2; CD325; CDHN; CDw325; NCAD ;

Immunogen:

Human N-cadherin.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A